



Product Datasheet

COX4-1 antibody (orb323124)

Description	Mouse monoclonal antibody to COX4-1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant protein corresponding to human COX4-1.
Target	COX4I1
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Mouse
Uniprot ID	P19783 , P10888 , P13073
Entrez	29445 , 12857 , 1327
Dilution Range	WB: 1:1000-1:3000, IF/ICC: 1:100-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of COX4-1 expressi...



Immunofluorescent analysis of COX4-1 sta...