



Product Datasheet

Cytokeratin 7 antibody (orb323133)

Description	Mouse monoclonal antibody to Cytokeratin 7
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IH, IP, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence of human Cytokeratin 7. The exact sequence is proprietary.
Target	KRT7
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Mouse
Uniprot ID	Q6IG12, Q9DCV7, P08729
Entrez	110310, 300242, 3855
Dilution Range	WB: 1:1000-1:2000, IHC-P: 1:100-1:200, IP: 1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of Cytokeratin 7 e...



Immunohistochemical analysis of Cytokera...



Immunoprecipitation of Cytokeratin 7 fro...