



Product Datasheet

EGFR antibody (orb323136)


Description	Mouse monoclonal antibody to EGFR
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence of human EGFR. The exact sequence is proprietary.
Target	EGFR
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Mouse
Uniprot ID	P00533
Entrez	1956
Dilution Range	WB: 1:1000-1:3000, IHC-P: 1:200-1:500, IF/ICC: 1:100-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of EGFR expression...



Immunohistochemical analysis of EGFR staining...



Immunofluorescent analysis of EGFR staining...