


Product Datasheet

GAPDH antibody (orb323147)

Description	Mouse monoclonal antibody to GAPDH
Species/Host	Mouse
Reactivity	Human, Insect, Mouse, Porcine, Primate, Rabbit, Rat, Sheep, Yeast
Conjugation	Unconjugated
Tested Applications	IH, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence of human GAPDH. The exact sequence is proprietary.
Target	GAPDH
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 0.2% rAlbumin, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 0.2% rAlbumin, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Mouse
Uniprot ID	P16858, P04406, P04797
Entrez	100042025, 2597, 24383
Dilution Range	WB: 1:2000-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of GAPDH expressio...