



Product Datasheet

HSC70 antibody (orb323171)

Description	Mouse monoclonal antibody to HSC70
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IH, WB
Immunogen	Recombinant protein corresponding to human HSC70.
Target	HSPA8
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Mouse
Uniprot ID	P11142
Entrez	3312
Dilution Range	WB: 1:2000-1:5000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of HSC70 expression...



Immunohistochemical analysis of HSC70 staining...