



Product Datasheet

HSP90 alpha antibody (orb323174)

Description	Mouse monoclonal antibody to HSP90 alpha
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IH, WB
Immunogen	Recombinant protein corresponding to human HSP90 alpha.
Target	HSP90AA1
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Mouse
Uniprot ID	P07900
Entrez	3320
Dilution Range	WB: 1:2000-1:5000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of HSP90 alpha exp...



Immunohistochemical analysis of HSP90 al...