



Product Datasheet

NF-kappaB p65 antibody (orb323181)


Description	Mouse monoclonal antibody to NF-kappaB p65
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IH, IP, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence of human NF-kappaB p65. The exact sequence is proprietary.
Target	RELA
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Mouse
Uniprot ID	Q04207 , Q04206
Entrez	19697 , 5970
Dilution Range	WB: 1:1000-1:3000, IHC-P: 1:100-1:200, IP: 1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of NF-kappaB p65 e...



Immunohistochemical analysis of NF-kappa...



Immunoprecipitation of NF-kappaB p65 fro...