



Product Datasheet

Peroxiredoxin 1 antibody (orb323187)

Description	Mouse monoclonal antibody to Peroxiredoxin 1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant protein corresponding to human Peroxiredoxin 1.
Target	PRDX1
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Mouse
Uniprot ID	P35700, Q63716, Q06830
Entrez	5052, 117254, 18477
Dilution Range	WB: 1:1000-1:3000, IF/ICC: 1:100-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of Peroxiredoxin 1...



Immunofluorescent analysis of Peroxiredo...