

Product Datasheet

NR2F2 antibody (orb32322)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	NR2F2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	45571
Uniprot ID	P24468
NCBI	NP_001138629.1 , NP_001138628.1 , NP_066285.1 , NP_001138627.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

