



## Product Datasheet

### c-FOS antibody (orb323248)

|                            |                                                                                                           |
|----------------------------|-----------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to c-FOS                                                                       |
| <b>Species/Host</b>        | Rabbit                                                                                                    |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                                                         |
| <b>Conjugation</b>         | Unconjugated                                                                                              |
| <b>Tested Applications</b> | IH, WB                                                                                                    |
| <b>Immunogen</b>           | Recombinant protein corresponding to human c-FOS.                                                         |
| <b>Target</b>              | FOS                                                                                                       |
| <b>Preservatives</b>       | Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide. |
| <b>Form/Appearance</b>     | Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide. |
| <b>Storage</b>             | Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.          |
| <b>Note</b>                | For research use only                                                                                     |
| <b>Clonality</b>           | Polyclonal                                                                                                |
| <b>Source</b>              | Rabbit                                                                                                    |
| <b>Uniprot ID</b>          | <a href="#">P01100</a> , <a href="#">P12841</a> , <a href="#">P01101</a>                                  |
| <b>Entrez</b>              | <a href="#">2353</a> , <a href="#">14281</a> , <a href="#">140675</a>                                     |
| <b>Dilution Range</b>      | WB: 1:1000-1:2000, IHC-P: 1:200-1:500                                                                     |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                           |



Western blot analysis of c-FOS expression...



Immunohistochemical analysis of c-FOS staining...