


## Product Datasheet

### His-tag antibody (HRP) (orb323283)

|                            |                                                                                                                     |
|----------------------------|---------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody against His-tag conjugated to HRP                                                         |
| <b>Species/Host</b>        | Rabbit                                                                                                              |
| <b>Conjugation</b>         | HRP                                                                                                                 |
| <b>Tested Applications</b> | WB                                                                                                                  |
| <b>Immunogen</b>           | KLH-conjugated synthetic peptide encompassing a sequence of His-tag-HRP labeled. The exact sequence is proprietary. |
| <b>Preservatives</b>       | Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.           |
| <b>Form/Appearance</b>     | Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.           |
| <b>Storage</b>             | Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.                    |
| <b>Note</b>                | For research use only                                                                                               |
| <b>Clonality</b>           | Monoclonal                                                                                                          |
| <b>Source</b>              | Mouse                                                                                                               |
| <b>Dilution Range</b>      | WB: 1:2000-1:5000                                                                                                   |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                     |



Western blot analysis of over-expressed ...