


Product Datasheet

Myc-tag antibody (HRP) (orb323284)

Description	Mouse monoclonal antibody against Myc-tag conjugated to HRP
Species/Host	Rabbit
Conjugation	HRP
Tested Applications	WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence of Myc-tag-HRP labeled. The exact sequence is proprietary.
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Mouse
Dilution Range	WB: 1:1000-1:3000
Expiration Date	12 months from date of receipt.



Western blot analysis of over-expressed ...