



Product Datasheet

Beta-tubulin antibody (HRP) (orb323286)

Description	Mouse monoclonal antibody against Beta-tubulin conjugated to HRP
Species/Host	Mouse
Reactivity	A. thaliana, Human, Mouse, Porcine, Primate, Rabbit, Rat, Zebrafish
Conjugation	HRP
Tested Applications	WB
Immunogen	Recombinant protein corresponding to human Beta-tubulin.
Target	TUBB
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Mouse
Uniprot ID	P07437
Entrez	203068
Dilution Range	WB: 1:20000-1:600000, IHC-P: 1:1000-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of Beta-tubulin An...



Western blot analysis of Beta-tubulin An...