

Product Datasheet

MCP1 antibody (orb323291)

Description

MCP1 antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Porcine, Rat

Conjugation

Unconjugated

Tested

ELISA, ICC, IF, IHC-P, WB

Applications

Immunogen

KLH conjugated synthetic peptide derived from human MCP1. Please contact us for the exact immunogen sequence. The peptide is available as orb12927.

Target

MCP1. Recognizes endogenous levels of MCP1 protein

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

11 kDa

Uniprot ID

[P13500](#), [P10148](#), [P14844](#)

NCBI

[002973](#), [002982](#)

Entrez


[6347](#)

Dilution Range


IHC-P: 1:100-800, WB: 1:100-500, IF/ICC: 1:100-800

Expiration Date


12 months from date of receipt.



IHC-P staining of human thyroid cancer t...



IHC-P image of pig large intestines tiss...



Immunohistochemical staining of rat brai...