



Product Datasheet

EpCAM antibody (orb323292)


Description	EpCAM antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human EpCAM. Please contact us for the exact immunogen sequence. The peptide is available as orb374948.
Target	EpCAM. Recognizes endogenous levels of EpCAM protein
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	35 kDa
Uniprot ID	O55159, Q99JW5, P16422
Entrez	4072
Dilution Range	WB: 1:100-300, IHC-P: 1:100-500 (based on 0.5 mg/ml)
Expiration Date	12 months from date of receipt.



WB analysis of rat liver (lane 1), human...



Western blot analysis of rat liver (lane...)



IHC-P staining of human endometrial cancer...