


Product Datasheet

PYHIN1 antibody (orb32390)

Description	PYHIN1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Pyrin and HIN domain-containing protein 1 protein (1-245AA)
Target	PYHIN1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q6K0P9
Dilution Range	WB:1:1000-5000
Expiration Date	12 months from date of receipt.



Western blot
analysis of
Mouse small
in...