



Product Datasheet

MIK1 antibody (orb32399)


Description	MIK1 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human MIK1. Please contact us for the exact immunogen sequence. The peptide is available as orb374793.
Target	MIK1
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	54 kDa
Uniprot ID	Q8NB16 , Q9D2Y4 , D3ZKP6
NCBI	001135969 , 001142497
Entrez	197259
Dilution Range	WB: 1:200-2000, IHC-P: 1:100-500, IF/ICC: 1:100-500
Expiration Date	12 months from date of receipt.



IHC-P staining of human breast tumour ti...



Immunohistochemical staining of rat skin...



IHC-P image of human thyroid tumour tiss...