



Product Datasheet

ARNT antibody (orb32444)

Description	Purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Synthetic Peptide
Target	ARNT
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	86636
Uniprot ID	P27540
NCBI	NP_848514.1 , NP_001659.1 , NP_001184254.1
Dilution Range	IF/ICC: 1:10-50, WB: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of NCI-H460 cell lysate.