


Product Datasheet

NKX2-3 antibody (orb324447)

Description	Rabbit polyclonal antibody to NKX2-3
Species/Host	Rabbit
Reactivity	Canine, Equine, Guinea pig, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the N terminal region of human NKX2-3
Target	NKX2-3
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	1.0 mg/ml
Storage	For short term use, store at 2-8C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: ADLEHHFHSAPCMLAAAEQTQFSDGGEEDEEDEGEKLSYLSLAAADGHG
Clonality	Polyclonal
MW	32kDa
Uniprot ID	Q8TAU0
NCBI	NP_660328
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse Spleen ti...



Western blot analysis of human Spleen ti...