


Product Datasheet

CLDN17 antibody (orb324486)


Description	Rabbit polyclonal antibody to CLDN17
Species/Host	Rabbit
Reactivity	Canine, Equine, Guinea pig, Human, Mouse, Porcine, Rat, Sheep
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of human CLDN17
Target	CLDN17
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	1.0 mg/ml
Storage	For short term use, store at 2-8C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: KQVQCTGSNERAKAYLLGTSGVLFILTGIFVLIPVSWTANIIIRDFYNPA
Clonality	Polyclonal
MW	25kDa
Uniprot ID	P56750
NCBI	NP_036263
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human Sp...



Western blot analysis of Jurkat cell lys...



Immunohistochemical staining of human Lu...