



Product Datasheet

AKAP17A antibody (orb324519)

Description	Rabbit polyclonal antibody to AKAP17A
Species/Host	Rabbit
Reactivity	Canine, Equine, Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	The immunogen is a synthetic peptide directed towards the N terminal region of human DXYS155E
Target	AKAP17A
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: NWEVMERLKGMVQNHQFSTLRISKSTMDFIRFEGEVENKSLVKSFLACLD
Clonality	Polyclonal
MW	81kDa
Uniprot ID	Q8N6U9
NCBI	NP_005079
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat cell lys...



Immunohistochemical staining of human ki...