



Product Datasheet

SUV420H1 antibody (orb324552)


Description	Rabbit polyclonal antibody to SUV420H1
Species/Host	Rabbit
Reactivity	Canine, Equine, Guinea pig, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	The immunogen is a synthetic peptide directed towards the C terminal region of human SUV420H1
Target	KMT5B
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	1.0 mg/ml
Storage	For short term use, store at 2-8C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: NNGFNSGSGSSSTKLKIQKRDEENRGSYTEGLHENGVCSDPLSLLES
Clonality	Polyclonal
MW	52kDa
Uniprot ID	Q4V775
NCBI	AAH98121
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human Sk...



Western blot analysis of human Placenta ...



Western blot analysis of Jurkat Whole Ce...