



Product Datasheet

SLC18A1 antibody (orb324613)


Description	Rabbit polyclonal antibody to SLC18A1
Species/Host	Rabbit
Reactivity	Equine, Guinea pig, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of human SLC18A1
Target	SLC18A1
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: YYLRSPPAKEEKLAILSQDCPMETRMATQKPTKEFPLGEDSDEEPPDHEE
Clonality	Polyclonal
MW	56 kDa
Uniprot ID	P54219
NCBI	NP_003044
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human brain tissue.



Western blot analysis of 293T cell lysate.



Immunohistochemical staining of human brain tissue.