


Product Datasheet

VEZF1 antibody (orb324621)


Description	Rabbit polyclonal antibody to VEZF1
Species/Host	Rabbit
Reactivity	Canine, Equine, Guinea pig, Human, Mouse, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of human VEZF1
Target	VEZF1
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8C up to 1 week. For long term storage, store at -20C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q14119
Expiration Date	12 months from date of receipt.

VEZF1




Western blot analysis of human Fetal Liv...

VEZF1



Western blot analysis of human Fetal Lun...

VEZF1



Western blot analysis of human Fetal Hea...