


Product Datasheet


PRDM9 antibody (orb324624)

Description	Rabbit polyclonal antibody to PRDM9
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of human PRDM9
Target	PRDM9
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: CPGDQNQEQQYPDPHSRNDKTKGQEIKERSKLLNKRTWQREISRAFSSPP
Clonality	Polyclonal
MW	103kDa
Uniprot ID	Q9NQV7
NCBI	NP_064612
Expiration Date	12 months from date of receipt.


PRDM9



Western blot analysis of human HepG2 tis...



Western blot analysis of 293T cells lysa...



Immunohistochemical staining of human Te...