

Product Datasheet

9330134C04Rik antibody [Out of stock] (orb324644)

Description Rabbit polyclonal antibody to 9330134C04Rik

Species/Host Rabbit

Reactivity Canine, Equine, Guinea pig, Human, Mouse, Porcine, Rat

Conjugation Unconjugated

Tested Applications WB

Target 9330134C04Rik

Preservatives Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance Liquid

Concentration 0.5 mg/ml

Storage For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.

Note For research use only

Clonality Polyclonal

Uniprot ID **Q6DI94**

Expiration Date 12 months from date of receipt.



Western
blot
analysis of
mouse
Liver tis...