


Product Datasheet

CACNB4 antibody (orb324660)

Description	Rabbit polyclonal antibody to CACNB4
Species/Host	Rabbit
Reactivity	Canine, Equine, Guinea pig, Human, Mouse, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of human CACNB4
Target	CACNB4
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: FDGRISITRVTDISLAKRSVLNPNPKRAIIERSNTRISSLAEVQSEIERIF
Clonality	Polyclonal
MW	50kDa
Uniprot ID	A7BJ74
NCBI	NP_001139270
Expiration Date	12 months from date of receipt.

CACNB2




Western blot analysis of human Placenta ...



Western blot analysis of MCF7 cell lysat...

CACNB2



Western blot analysis of human Placenta ...