

Product Datasheet

Cphx1 antibody (orb324718)

Description	Rabbit polyclonal antibody to Cphx1
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the N terminal region of mouse C330003B14RIK
Target	Cphx1
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	1.0 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: MLSKNFPGAPETKDNRSKARKRYGSRNSKPRHKFSRDELKRLKQEFAYAP
Clonality	Polyclonal
MW	20kDa
Uniprot ID	Q8BX39
NCBI	NP_780551
Expiration Date	12 months from date of receipt.



Western blot analysis of SP2/0 cell lysa...