

Product Datasheet

NR1H5 antibody [Out of stock] (orb324728)

Description	Rabbit polyclonal antibody to NR1H5
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the N terminal region of mouse NR1H5
Target	NR1H5
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	1.0 mg/ml
Storage	For short term use, store at 2-8C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q811W9
Expiration Date	12 months from date of receipt.



Western blot analysis of SP2/0 cell lysis...