


Product Datasheet

ZSCAN29 antibody (orb324741)

Description	Rabbit polyclonal antibody to ZSCAN29
Species/Host	Rabbit
Reactivity	Canine, Equine, Guinea pig, Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen for Anti-ZSCAN29 antibody is: synthetic peptide directed towards the middle region of Human ZSC29
Target	ZSCAN29
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: GRPRSSVTSVKQGQEVRLKMTTPPKSSQELLSVRQESVEPQPRGVPPKKER
Clonality	Polyclonal
MW	50 kDa
Uniprot ID	Q8IWIY8
Expiration Date	12 months from date of receipt.



Western blot analysis of human RPMI-8226...