


Product Datasheet

ZNF419 antibody (orb324808)

Description	Rabbit polyclonal antibody to ZNF419
Species/Host	Rabbit
Reactivity	Canine, Equine, Human, Mouse, Porcine, Rat, Sheep
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the C terminal region of human ZNF419A
Target	ZNF419
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	1.0 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: GRLFRENSSLVKHQRVHTGAKPYECRECGKFFRHNSSLFKHRRHTGEMQ
Clonality	Polyclonal
MW	59kDa
Uniprot ID	Q96HQ0
NCBI	NP_078967
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2 cell lysa...