



Product Datasheet

NACC1 antibody (orb324816)


Description	Rabbit polyclonal antibody to NACC1
Species/Host	Rabbit
Reactivity	Canine, Equine, Guinea pig, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	The immunogen is a synthetic peptide directed towards the C terminal region of human BTBD14B
Target	NACC1
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	1.0 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: WMPKVKVLKAEDDAYTTFISETGKIEPDMMGVHGFETASHEGEAGPSAE
Clonality	Polyclonal
MW	57kDa
Uniprot ID	Q96RE7
NCBI	NP_443108
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human kidney tissue.



Western blot analysis of HepG2 cell lysate.



Immunohistochemical staining of human liver tissue.