

Product Datasheet

MUPCDH antibody (orb32482)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	CDHR5
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	88223
Uniprot ID	Q9HBB8
NCBI	NP_068743.2 , NP_112554.2 , NP_001165439.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

