

Product Datasheet

ZNF709 antibody (orb324825)

Description	Rabbit polyclonal antibody to ZNF709
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the N terminal region of human ZNF709
Target	ZNF709
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	1.0 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: DVMQETFVNLASIGENWEEKNIEDHKNQGRKLRSHMVERLCERKEGSQFG
Clonality	Polyclonal
MW	75kDa
Uniprot ID	Q8N972
NCBI	NP_689814
Expiration Date	12 months from date of receipt.

