

Product Datasheet

GRSF1 antibody (orb324854)


Description	Rabbit polyclonal antibody to GRSF1
Species/Host	Rabbit
Reactivity	Canine, Equine, Guinea pig, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of human GRSF1
Target	GRSF1
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: IRNGENGIHFLNDRDGKRRGDALIEMESEQDVQKALEKHRMYMGQRYVEV
Clonality	Polyclonal
MW	53kDa
Uniprot ID	Q12849
NCBI	NP_002083
Expiration Date	12 months from date of receipt.

GRSF1



Western blot analysis of human Jurkat ti...

GRSF1



Western blot analysis of human Hela tiss...

GRSF1



Western blot analysis of human Fetal Bra...