



Product Datasheet

TMEM63A antibody (orb324898)

Description	Rabbit polyclonal antibody to TMEM63A
Species/Host	Rabbit
Reactivity	Guinea pig, Human, Mouse, Rat, Yeast
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the N terminal region of human TMEM63A
Target	TMEM63A
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: MDSPFLELWQSKAVSIREQLGLGDRPNDSYCYNSAKNSTVLQGVTFGGIP
Clonality	Polyclonal
MW	89kDa
Uniprot ID	O94886
NCBI	NP_055513
Expiration Date	12 months from date of receipt.



Western blot analysis of 293T cell lysat...



Western blot analysis of human stomach(T...