



Product Datasheet

GTF2B antibody (orb32495)

Description	Rabbit polyclonal antibody to GTF2B (General transcription factor IIB) is also known as TFIIIB, is on...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	GTF2B
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	34833
Uniprot ID	Q00403
NCBI	NP_001505.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse kidney ti...