

Product Datasheet

WRAP73 antibody (orb324950)

Description	Rabbit polyclonal antibody to WRAP73
Species/Host	Rabbit
Reactivity	Equine, Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of human WDR8
Target	WRAP73
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8C up to 1 week. For long term storage, store at -20C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q9P2S5
Expiration Date	12 months from date of receipt.



Western blot analysis of human Thymus ti...