



Product Datasheet

WDR4 antibody (orb324953)

Description	Rabbit polyclonal antibody to WDR4
Species/Host	Rabbit
Reactivity	Equine, Human, Mouse, Porcine, Rat, Yeast
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	The immunogen is a synthetic peptide directed towards the C terminal region of human WDR4
Target	WDR4
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	1.0 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: AGADASFSSLYKATFDNVTSYLKKKEERLQQLEKKQRRRSPPPGPDGHA
Clonality	Polyclonal
MW	45kDa
Uniprot ID	P57081
NCBI	NP_387510
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2 cell lysa...



Immunohistochemical staining of human Ki...