



Product Datasheet

KRT2 antibody (orb324963)


Description	Rabbit polyclonal antibody to KRT2
Species/Host	Rabbit
Reactivity	Canine, Equine, Guinea pig, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of human KRT2A
Target	KRT2
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8C up to 1 week. For long term storage, store at -20C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P35908
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human Lu...



Western blot analysis of Jurkat cell lys...



Immunohistochemical staining of human In...