



Product Datasheet

SLC7A7 antibody (orb325097)


Description	Rabbit polyclonal antibody to SLC7A7
Species/Host	Rabbit
Reactivity	Canine, Equine, Guinea pig, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of human SLC7A7
Target	SLC7A7
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: WGTLVQDIFTYAKVLALIGIVRLGQGASTHFENSFEGSSFAVGDI
Clonality	Polyclonal
MW	56kDa
Uniprot ID	Q9UM01
NCBI	NP_003973
Expiration Date	12 months from date of receipt.



Western blot analysis of PANC1 cell lysa...



Western blot analysis of human Fetal Hea...



Western blot analysis of human Fetal Lun...