

Product Datasheet

SLC7A8 antibody (orb325105)

Description	Rabbit polyclonal antibody to SLC7A8
Species/Host	Rabbit
Reactivity	Canine, Equine, Guinea pig, Human, Mouse, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of human SLC7A8
Target	SLC7A8
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	1.0 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: PRAIFISIPLVTFVYVFANVAYVTAMSPQELLASNAVAVTFGEKLLGVMA
Clonality	Polyclonal
MW	37kDa
Uniprot ID	Q9UHI5
NCBI	NP_877392
Expiration Date	12 months from date of receipt.



Western blot analysis of Fetal Liver tis...