



Product Datasheet

SLC46A3 antibody (orb325132)

Description	Rabbit polyclonal antibody to SLC46A3
Species/Host	Rabbit
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	The immunogen is a synthetic peptide directed towards the N terminal region of human SLC46A3
Target	SLC46A3
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	1.0 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: MKILFVEPAIFLSAFAMTLTGPLTTQYVYRRIWEETGNYTFSSDSNISEC
Clonality	Polyclonal
MW	51kDa
Uniprot ID	Q7Z3Q1
NCBI	NP_861450
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat cell lysates



Immunohistochemical staining of human kidney tissue