


Product Datasheet

TSEAR antibody (orb325240)

Description	Rabbit polyclonal antibody to TSEAR
Species/Host	Rabbit
Reactivity	Canine, Equine, Guinea pig, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the C-terminal region of Human TSEAR
Target	TSPEAR
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: YDWEFFSVGPYSFLVVANTFNGTSTKVHSHLYIRLLGSFQLFQSFTFGA
Clonality	Polyclonal
MW	73kDa
Uniprot ID	Q8WU66
NCBI	NP_659428
Expiration Date	12 months from date of receipt.



Western blot analysis of human 293T cells transfected with TSEAR. The blot shows a single band at approximately 73 kDa, indicating successful expression of the target protein.