

Product Datasheet

GPR137 antibody (orb32532)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	GPR137
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	46059
Uniprot ID	Q96N19
NCBI	NP_001164197.1, NP_001170829.1, NP_001164351.1, NP_064540.3, NP_001164352.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

