


Product Datasheet

exd antibody (orb325377)


Description	Rabbit polyclonal antibody to exd
Species/Host	Rabbit
Reactivity	Canine, Equine, Guinea pig, Human, Mouse, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide corresponding to a region of Fruit fly
Target	exd
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	1.0 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: EAKEELARKCGITVSQVSNWFGNKRIRYKKNIGKAQEEANLYAAKKAAGA
Clonality	Polyclonal
MW	42kDa
Uniprot ID	P40427
NCBI	NP_727924
Expiration Date	12 months from date of receipt.



Western blot analysis of exd antibody



Western blot analysis of exd antibody



Western blot analysis of r-Exd tissue us...