


Product Datasheet

TTC7A antibody (orb325425)

Description	Rabbit polyclonal antibody to TTC7A
Species/Host	Rabbit
Reactivity	Canine, Equine, Guinea pig, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the N-terminal region of human TTC7A
Target	TTC7A
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: EAGEFLPKGYLALGLTYSLQATDATLKSQDELHRKALQTLERAQQLAPS
Clonality	Polyclonal
MW	55kDa
Uniprot ID	Q9ULT0
NCBI	NP_001275880.1
Expiration Date	12 months from date of receipt.



Western blot analysis of human Stomach T...