


Product Datasheet

CGAT2 antibody (orb325473)

Description	Rabbit polyclonal antibody to CGAT2
Species/Host	Rabbit
Reactivity	Canine, Equine, Guinea pig, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the N-terminal region of Human CGAT2
Target	CSGALNACT2
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: PGVVGENYGKEYYQALLQEQQEEHYQTRATSLKRQIAQLKQELQEMSEKMR
Clonality	Polyclonal
MW	59kDa
Uniprot ID	Q8N6G5
NCBI	XP_005271877
Expiration Date	12 months from date of receipt.



Western blot analysis of human HepG2 cells. The band at approximately 59 kDa represents CGAT2 protein.