


Product Datasheet

UNC93B1 antibody (orb325532)

Description	Rabbit polyclonal antibody to UNC93B1
Species/Host	Rabbit
Reactivity	Canine, Equine, Guinea pig, Human, Mouse, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the N terminal region of human UNC93B1
Target	UNC93B1
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	1.0 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: LKNVLAASAGGMLTYGVYLGLLQMQLILHYDETYREVKYGNMGLPDIDSK
Clonality	Polyclonal
MW	67kDa
Uniprot ID	Q9H1C4
NCBI	NP_112192
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat cell lys...