


Product Datasheet

LN28B antibody (orb325570)

Description	Rabbit polyclonal antibody to LN28B
Species/Host	Rabbit
Reactivity	Canine, Equine, Guinea pig, Human, Mouse, Porcine, Rat, Sheep
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the C-terminal region of Human LN28B
Target	LIN28B
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: GCTSPFPQEAREISERSGRSPQEASSTKSSIAPEEQSKKGPSVQKRKK
Clonality	Polyclonal
MW	27kDa
Uniprot ID	Q6ZN17
NCBI	NP_001004317
Expiration Date	12 months from date of receipt.



Western blot analysis of human 786-0 cells. The blot shows a single band at approximately 27 kDa, indicating the presence of the target protein LN28B.